



# SARS-CoV-2 Antigen Rapid Test ( Saliva )

**Speicheltest**

## INTENDED USE

SARS-CoV-2 Antigen Rapid Test (Saliva) is intended for the qualitative detection of SARS-CoV-2 Antigen in human saliva specimens in vitro.


## FEATURES

- Easy to operate, no other equipment needed.
- High Accuracy, Specificity and Sensitivity.
- Show results within 20 minutes.
- Room temperature preservation.




[www.wizbiotech.com](http://www.wizbiotech.com)


## OPERATION PROCESS




1. Gently spit saliva into the disposable sample container until the liquid saliva (non-bubble) reaches about 0.5cm depth of the container.




2. Hold a disposable sampling rod, and use its swab head to stir the saliva counter clockwise / clockwise (about 5 circles).




3. Loosen the lower cover (the extraction tube is composed of tube, lower cover and upper cover).




4. Immerse the swab head into the extraction solution in the extraction tube, and rotate it close to the wall of the sample extraction tube for about 10 seconds or 10 times to make the sample dissolve in the solution as much as possible.



5. Squeeze the swab head along the inner wall of the sample extraction tube to keep the liquid in the tube as much as possible, discard the disposable sampling rod or break the disposable sampling rod to leave the swab head in the extraction tube, tighten the cover, gently shake the extraction tube, gently mix and stand by.



6. Open the upper cover of the extraction tube (the extraction tube with processed specimens), add 3 drops vertically into the sample well of the test device.




7. The test results should be interpreted within 20 minutes, the result is invalid if more than 30 minutes.

## LIMIT OF DETECTION

$1.7 \times 10^2$  TCID<sub>50</sub>/mL


## INTERPRETATION OF RESULT



Positive



Negative



Invalid

## CLINICAL PERFORMANCE

(Saliva samples)

Test Results of Saliva samples	Reference PCR Results		
	POS	NEG	Total
POS	328	0	328
NEG	16	555	571
Total	344	555	899

PPA: 95.35% (C.I. 92.58%~97.12%) NPA: 100.00% (C.I. 99.31%~100.00%) OPA: 98.22% (C.I. 97.13%~98.90%)

## SPECIFICATION

25Tests/Kit, 1Test/Kit



Xiamen Wiz Biotech Co., Ltd. Web: [www.wizbiotech.com](http://www.wizbiotech.com)

- Tel: +86-592-6808278 Fax: +86-592-6808279 E-mail: [sales@wizbiotech.com](mailto:sales@wizbiotech.com)
- Address: 3-4 Floor, NO.16 Building, Bio-medical Workshop, 2030 Wengjiao Xi Road, Haicang District, Xiamen City, Fujian Province, 361026, P.R. China



## EC Declaration of Conformity

**Manufacturer** Xiamen Wiz Biotech Co., Ltd.  
3-4 Floor, NO.16 Building, Bio-medical Workshop, 2030 Wengjiao Xi Road, Haicang District, Xiamen City, Fujian Province, 361026, P.R. China

**European Representative** Qarad EC-REP BV  
Pas 257, 2440 Geel, Belgium

**Product** SARS-CoV-2 Antigen Rapid Test (Saliva)

**Model** 1 test/kit, 2 tests/kit, 20 tests/kit, 25 tests/kit

**Catalogue number** 51232701, 51232702, 51232703, 51232704

**Classification** Others

**Conformity assessment route:** Annex III (IVDD 98/79 EC)


We, the manufacturer, herewith, declares that the product(s) as specified above meet(s) the applicable provisions of the European Directive 98/79/EC on *in vitro* Diagnostic Medical Devices. All supporting technical documentation demonstrating compliance is kept by the manufacturer and can be made available by the Authorized Representative in Europe.

**General applicable directive:**

Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on *in vitro* diagnostic medical devices.

Signed on 04/(Day) 02/(Month) of 2021. Place Xiamen.

**Represented by**

Signature:   
Name of authorized signatory: Xiaozhen Wang  
Position held in the company: General Manager

Seal/Stamp:



Xiamen Wiz Biotech Co., Ltd. Web: [www.wizbiotech.com](http://www.wizbiotech.com)

- Tel: +86-592-6808278 Fax: +86-592-6808279 E-mail: [sales@wizbiotech.com](mailto:sales@wizbiotech.com)
- Address: 3-4 Floor, NO.16 Building, Bio-medical Workshop, 2030 Wengjiao Xi Road, Haicang District, Xiamen City, Fujian Province, 361026, P.R. China



## SARS-CoV-2 Antigen Rapid Test (Sputum/Saliva/Stool)

### Clinical Report

The performance of the SARS-CoV-2 Antigen Rapid Test (Sputum/Saliva/Stool) was established with collected from 1012 individual symptomatic patients who were suspected of COVID-19. It contains 1012 sputum samples, 899 saliva samples and 899 stool samples. The collection, treatment, storage, transportation and detection of samples met the relevant requirements of the Instruction for Use. At the same time, SARS-CoV-2 was detected by emergency nucleic acid detection reagent.

Table 1 .Performance summary of the WIZ's SARS-CoV-2 Antigen Rapid Test(Sputum/Saliva/Stool)  
(Sputum samples)

WIZ Results	Reference PCR Results		
	POS	NEG	Total
POS	407	0	407
NEG	14	591	605
Total	421	591	1012

PPA: 96.67% (C.I. 94.50%~98.01%)

NPA: 100.00% (C.I. 99.35%~100.00%)

OPA: 98.62% (C.I. 97.69%~99.17%)

Table 2 .Performance summary of the WIZ's SARS-CoV-2 Antigen Rapid Test(Sputum/Saliva/Stool)  
(Saliva samples)

Test Results of Saliva samples	Reference PCR Results		
	POS	NEG	Total
POS	328	0	328
NEG	16	555	571
Total	344	555	899

PPA: 95.35% (C.I. 92.58%~97.12%)

NPA: 100.00% (C.I. 99.31%~100.00%)

OPA: 98.22% (C.I. 97.13%~98.90%)

Table 3. Performance summary of the WIZ's SARS-CoV-2 Antigen RapidTest(Sputum/Saliva/Stool)

(Stool samples)

Test Results of Stool samples	Reference PCR Results		
	POS	NEG	Total
POS	330	0	330
NEG	14	555	569
Total	344	555	899

PPA: 95.93% (C.I. 93.29%~97.56%)

NPA: 100.00% (C.I. 99.31%~100.00%)

OPA: 98.44% (C.I. 97.40%~99.07%)

## EXPLANATION OF TERMS:


C.I.: Confidence Interval

PPA: Positive Percent Agreement=True Positives / (True Positives +False Negatives)

NPA: Negative Percent Agreement=True Negatives / (True Negatives +False Positive)

OPA: Overall Percent Agreement= (True Positives +True Negatives) /Total Samples







CE

# SARS-CoV-2 Antigen Rapid Test (Saliva)

Speicheltest

Awarding ISO13485 Certification



REF  
LOT  
M

20°C 30°C



QC PASS  
WIZ BIOTECH

## SARS-CoV-2 Antigen Rapid Test (Saliva) Speicheltest

CE





## EC Declaration of Conformity

**Manufacturer** Xiamen Wiz Biotech Co., Ltd.  
3-4 Floor, NO.16 Building, Bio-medical Workshop, 2030 Wengjiao Xi  
Road, Haicang District, Xiamen City, Fujian Province, 361026, P.R.  
China

**European Representative** Qarad EC-REP BV  
Pas 257, 2440 Geel, Belgium

**Product** SARS-CoV-2 Antigen Rapid Test (Saliva)

**Model** 1 test/kit, 2 tests/kit, 20 tests/kit, 25 tests/kit

**Catalogue number** 51232701, 51232702, 51232703, 51232704

**Classification** Others

**Conformity assessment route:** Annex III (IVDD 98/79 EC)


We, the manufacturer, herewith, declares that the product(s) as specified above meet(s) the applicable provisions of the European Directive 98/79/EC on *in vitro* Diagnostic Medical Devices. All supporting technical documentation demonstrating compliance is kept by the manufacturer and can be made available by the Authorized Representative in Europe.

### General applicable directive:

Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on *in vitro* diagnostic medical devices.

Signed on 04 /(Day) 02 /(Month) of 2021 . Place Xiamen .


### Represented by

Signature:   
Name of authorized signatory: Xiaozhen Wang  
Position held in the company: General Manager

Seal/Stamp:







CE

# SARS-CoV-2 Antigen Rapid Test

( Saliva )

*Speicheltest*

REF

LOT




Awarding ISO13485 Certification



QC PASS

WIZ BIOTECH



CE

# SARS-CoV-2 Antigen Rapid Test

( Saliva )

*Speicheltest*